
INSTRUCTION FOR USE (LATEX POWDER FREE)

STATEMENT AND CAUTION

- 1) This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals.
- 2) The penetration resistance and chemical resistance has **been assessed under laboratory conditions and relates only to the tested specimen** taken from the palm only (except in cases where the glove is equal to or over 400mm - where the cuff is tested also) and **relates only to the chemical tested**. It can be different if the chemical is used in a mixture.
- 3) It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- 4) When used, protective gloves may provide less resistance to the dangerous chemical due to changes in the physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves.
- 5) Before usage, inspect the gloves for any defects or imperfections.
- 6) There are no possible allergen within the glove that are known to cause harm to the wearer.

EU Type Examination and ongoing Conformity by Notified Body:- **CE 2777**

SATRA Technology Europe Ltd
Bracetown Business Park,
Clonee, D15 YN2P, Ireland.

EN ISO 374-1:2016 Classification Of Permeation Performance Level

| Measured Breakthrough time (min) | Permeation Performance Level |
|----------------------------------|------------------------------|
| > 10 | 1 |
| > 30 | 2 |
| > 60 | 3 |
| > 120 | 4 |
| > 240 | 5 |
| > 480 | 6 |

| Hand Size | 6 (XS) | 7 (S) | 8 (M) | 9 (L) | 10 (XL) |
|------------|--------|-------|-------|-------|---------|
| Min Length | 240mm | | | | |

Resistance against Bacteria and Fungi – **PASS**
Resistance against Virus - **PASS**

| Chemical Permeation (EN ISO 374-1:2016) | Level | Mean Degradation (EN 374-4:2013) | |
|---|-------|----------------------------------|--|
| K 40% Sodium Hydroxide | 3 | -15.8% | Degradation levels indicate the change in Puncture Resistance of the glove after exposure to the challenge chemical. |
| T 37% Formaldehyde | 3 | 42.5% | |
| P 30% Hydrogen Peroxide | 2 | 0.0% | |

EN 420:2003+A1:2009
Fit For Special Purpose
(Palm protection)

| | |
|----------------------------|---|
| CE 2777 UK CA MD | Comply to MDR 2017/745 Class I |
| | Comply to PPE REGULATION 2016/425 Cat III |
| | EN 455 – 1, 2, 3, 4 |

A powder-free patient examination glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner.

Caution:

This product contains natural rubber latex and other components which may cause allergic reaction including anaphylactic responses. Isolated case of such reactions to natural rubber latex or powder have been reported. After use, wearers should visually check the gloves and remove any contamination from the outer surface before removing the gloves from the hand. Alternatively, carefully peel the gloves off the hand so that the contaminated outer glove does not touch your skin.

Warning: Do not use this product if you have known allergy to chemical additives.

If you experience an allergic reactions to this product, discontinue use immediately and consult your physician.

Storage:

Store in a cool, dry place. In storage avoid excessive heat (40°C, 104°F). Opened box should be shielded from exposure to direct sun or fluorescent lighting.

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